



Dental Equipment
Manufacturer Ltd



Supporting Excellence in Dentistry

DemDetect™


Submerged Implant Detector



www.demdental.com

DemDetect™ Submerged Implant Detector Kit

- Fast, reliable and non-invasive
- Reduces risk of infection
- Faster treatment and recovery time
- Ergonomic design for easy operation
- Emits audio and visual signals to indicate implant location
- Autoclavable sensor heads
- Suitable for locating submerged implants in upper or lower jaw
- Builds practitioner confidence



EN

Kit Contains:

- 1** DemDetect™ Implant Detector - 1 pc
- 2** Sensor - 1 pc
- 3** Sensor holder - 1 pc
- 4** AA alkaline battery - 1 pc

SP

Contenido del Kit:

- 1** DemDetect™ Detector de Implante - 1 pc
- 2** Sensor - 1 pc
- 3** Soporte del sensor - 1 pc
- 4** AA batería alcalina - 1 pc

FR

Contenu du kit:

- 1** DemDetect™ Décteur d'Implant - 1 pc
- 2** Capteur - 1 pc
- 3** Support capteur - 1 pc
- 4** AA pile alkaline - 1 pc

DEM Dental Equipment Manufacturer Ltd manufactures and exports a range of dental products worldwide. We supply high-quality products for fast, safe, cost-effective dental treatment.

Quality Certification: ISO 13485



DEM Dental Equipment Manufacturer Ltd
Michael House, Castle Street, Exeter, EX4 3LQ, UK

Contact DEM's dedicated customer support team to discuss supply
dem@demdental.com | +44 (0)7543 213432

Distributed by: