






Dental Equipment
Manufacturer Ltd




Order Sample
Cat Number **FG-TC 7015L**




Dental Equipment
Manufacturer Ltd




Order Sample
Cat Number **FG-TC 7015L**



Dental Equipment
Manufacturer Ltd



Order Sample
Cat Number **FG-TC 7015L**



Carbide Operative	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
-------------------	------	---------------	-------------------	-----------	-----------

BALL

	TC 1/4	005	0.1	19.0	FG-TC1/4
	TC 1/2	006	0.5	19.0	FG-TC1/2
	TC 1	008	0.7	19.0	FG-TC1
	TC 2	010	0.8	19.0	FG-TC2
	TC 2SH short	010	0.8	16.0	FG-TC2SH
	TC 3	012	0.9	19.0	FG-TC3
	TC 4	014	1.1	19.0	FG-TC4
	TC 4SH short	014	1.1	16.0	FG-TC4SH
	TC 5	016	1.4	19.0	FG-TC5
	TC 6	018	1.6	19.0	FG-TC6
	TC 7	021	1.7	19.0	FG-TC7
	TC 8	023	1.9	19.0	FG-TC8

INVERTED CONE

	TC 331/2	006	0.5	19.0	FG-TC331/2
	TC 34	008	0.6	19.0	FG-TC34
	TC 35	010	0.8	19.0	FG-TC35
	TC 35SH short	010	0.8	16.0	FG-TC35SH
	TC 36	012	1.0	19.0	FG-TC36
	TC 37	014	1.1	19.0	FG-TC37

FLAT FISSURE Plain Cut

	TC 56	008	3.2	19.0	FG-TC56
	TC 56SH short	008	3.2	16.0	FG-TC56SH
	TC 57	010	3.8	19.0	FG-TC57
	TC 57SH short	010	3.8	16.0	FG-TC57SH
	TC 57L L-head	010	4.9	19.0	FG-TC57L
	TC 58	012	3.8	19.0	FG-TC58

Carbide Operative	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
-------------------	------	---------------	-------------------	-----------	-----------

TAPERED FISSURE Plain Cut

	TC 169	009	3.4	19.0	FG-TC169
	TC 169L L-head	009	4.7	19.0	FG-TC169L
	TC 170	010	3.8	19.0	FG-TC170
	TC 170L L-head	010	4.9	19.0	FG-TC170L

AMALGAM PREP

	TC 245	008	2.7	19.0	FG-TC245
	TC 245SH short	008	2.7	16.0	FG-TC245SH

PEAR

	TC 329	006	0.9	19.0	FG-TC329
	TC 330	008	1.9	19.0	FG-TC330
	TC 330S	008	2.0	19.0	FG-TC330S
	TC 330SH short	008	1.9	16.0	FG-TC330SH
	TC 331	010	1.8	19.0	FG-TC331
	TC 331SH short	010	1.8	16.0	FG-TC331SH
	TC 331L L-head	010	3.7	19.0	FG-TC331L
	TC 332	012	2.0	19.0	FG-TC332

FLAT FISSURE Cross Cut

	TC 556	008	3.2	19.0	FG-TC556
	TC 556SH short	008	3.2	16.0	FG-TC556SH
	TC 557	010	3.8	19.0	FG-TC557
	TC 557L	010	4.9	19.0	FG-TC557L
	TC 557SH short	010	3.8	16.0	FG-TC557SH
	TC 558	012	3.8	19.0	FG-TC558

Carbide Operative	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
-------------------	------	---------------	-------------------	-----------	-----------

TAPERED FISSURE Cross Cut

	TC 699	009	3.4	19.0	FG-TC699
	TC 700	010	3.8	19.0	FG-TC700
	TC 701	012	3.8	19.0	FG-TC701
	TC 701L L-head	012	4.9	19.0	FG-TC701L
	TC 702	016	4.3	19.0	FG-TC702

DOMED FISSURE Plain Cut

	TC 1156	080	3.2	19.0	FG-TC1156
	TC 1157	010	3.8	19.0	FG-TC1157
	TC 1157SH short	010	3.8	16.0	FG-T1157SH
	TC 1158	012	3.8	19.0	FG-TC1158

TAPERED FISSURE 8 Flutes

	TC 1171	012	3.8	19.0	FG-TC1171
--	---------	-----	-----	------	-----------

DOMED FISSURE Cross Cut

	TC 1556	008	3.2	19.0	FG-TC1556
	TC 1557	010	3.8	19.0	FG-C1557
	TC 1557SH short	010	3.8	16.0	FG-TC1557SH
	TC 1558	012	3.8	19.0	FG-TC1558

CROWN PREP

	TC 37923	023	4.3	19.0	FG-TC37923
	TC 84716	016	8.0	19.0	FG-TC84716
	TC 85616	016	8.0	19.0	FG-TC85616
	TC 85618	018	8.0	19.0	FG-TC85618

METAL/CROWN CUTTERS

	KC1	010	1.9	19.0	FG-KC1
	KC2	012	4.0	19.0	FG-KC2

Carbide Operative	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
-------------------	------	---------------	-------------------	-----------	-----------

METAL/CROWN CUTTERS

	TC 1557T	010	3.9	19.0	FG-TC 1557T
	TC 330T	008	2.0	19.0	FG-TC 330T
	TC 557T	010	3.7	19.0	FG-TC 557T

AESTHETIC FINISHERS 10 Straight Blades

	TC EF3	008	3.0	19.0	FG-TCEF3
	TC EF4	010	4.2	19.0	FG-TCEF4
	TC EF6	014	6.0	19.0	FG-TCEF6
	TC EF9	014	9.0	19.0	FG-TCEF9

AESTHETIC FINISHERS 12 Blades

	TC OS1	023	4.3	19.0	FG-TCOS1
	TC OS2	090	3.2	19.0	FG-TCOS2
	TC OS3	070	3.2	19.0	FG-TCOS3

TRIMMING & FINISHING 12 Blades

	TC 7006	018	1.6	19.0	FG-TC7006
	TC 7104	014	3.4	19.0	FG-TC7104
	TC 7106	018	3.7	19.0	FG-TC7106
	TC 274H	018	3.5	19.0	FG-TC274H
	TC 48L10H	010	8.0	19.0	FG-TC48L10H
	TC 48L12H	012	8.0	19.0	FG-TC48L12H
	TC 7404	014	3.3	19.0	FG-TC7404
	TC 7406	018	3.4	19.0	FG-TC 7406
	TC 7408	032	3.9	19.0	FG-TC7408
	TC 7642	010	7.0	19.0	FG-TC7642
	TC 7664	015	8.2	19.0	FG-TC7664
	TC 7675	016	8.0	19.0	FG-TC7675
	TC 7714	016	8.6	19.0	FG-TC7714
	TC 7802	010	3.4	19.0	FG-TC7802

Carbide Operative	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
-------------------	------	---------------	-------------------	-----------	-----------

TRIMMING & FINISHING 12 Blades

	TC 7803	012	3.4	19.0	FG-TC7803
	TC 7901	009	3.6	19.0	FG-TC7901
	TC 7902	010	3.6	19.0	FG-TC7902
	TC 7903	012	3.6	19.0	FG-TC7903

FINE TRIMMING & FINISHING 30 Blades

	TC EF9UF	014	9.0	19.0	FG-TCEF9UF
	TC 9103UF	012	3.6	19.0	FG-TC9103UF
	TC 9214UF	010	5.6	19.0	FG-TC9214UF
	TC 9406UF	018	3.4	19.0	FG-TC9406UF
	TC 9903UF	012	5.5	19.0	FG-TC9903UF

Carbide Surgical	Item	Head Ø1/10 mm	Head Length L1 mm	OAL mm	Order No.
------------------	------	---------------	-------------------	--------	-----------

BALL

	TC 1/2SL	060	0.6	25.0	FG-TC1/2SL
	TC 1SL	080	0.8	25.0	FG-TC1SL
	TC 2SL	010	1.0	25.0	FG-TC2SL
	TC 4SL	014	1.4	25.0	FG-TC4SL
	TC 6SL	018	1.8	25.0	FG-TC6SL
	TC 8SL	023	2.3	25.0	FG-TC8SL

LINDEMANN

	TC 161 Straight	010	1.6	25.0	FG-TC161
	TC 162 Cross Cut	010	1.6	25.0	FG-TC162

Carbide Surgical	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
------------------	------	---------------	-------------------	-----------	-----------

STRAIGHT FISSURE Cross Cut

	TC 557SL	010	3.8	25.0	FG-TC557SL
	TC 558SL	012	3.8	25.0	FG-TC558SL

TAPERED FISSURE Cross Cut

	TC 701SL	012	3.8	25.0	FG-TC701SL
	TC 702SL	016	4.3	25.0	FG-TC702SL


ENDO ACCESS

	151-Z	016	11.0	23.0	FG-151-Z
	151L-Z	016	11.0	28.0	FG-151L-Z
	152-016	016	9.0	23.0	FG-152-016






Dental Equipment
Manufacturer Ltd



Order Sample
Cat Number RA.TC0001 / HP.TC02SL



RA Carbide	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
------------	------	------------------	-------------------------	--------------	--------------

RA BALL / OPERATIVE

	RA.TC 1/2	006	0.5	22.0	RA.TC01/2
	RA.TC 1	008	0.7	22.0	RA.TC0001
	RA.TC 2	010	0.8	22.0	RA.TC0002
	RA.TC 3	012	0.9	22.0	RA.TC0003
	RA.TC 4	014	1.1	22.0	RA.TC0004
	RA.TC 5	016	1.4	22.0	RA.TC0005
	RA.TC 6	018	1.6	22.0	RA.TC0006
	RA.TC 7	021	1.7	22.0	RA.TC0007
	RA.TC 8	023	1.9	22.0	RA.TC0008

RA INVERTED CONE / OPERATIVE

	RA.TC 34	008	0.6	22.0	RA.TC0034
	RA.TC 35	010	0.8	22.0	RA.TC0035
	RA.TC 36	012	1.0	22.0	RA.TC0036

RA Carbide Surgical	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
------------------------	------	------------------	-------------------------	--------------	--------------

RA BALL / SURGICAL LENGTH

	RA.TC 2SL	010	1.0	26.0	RA.TC02SL
	RA.TC 4SL	014	1.4	26.0	RA.TC04SL
	RA.TC 6SL	018	1.8	26.0	RA.TC06SL
	RA.TC 8SL	023	2.3	26.0	RA.TC08SL

HP Carbide	Item	Head Ø1/10 mm	Head Length L1 mm	OAL L2 mm	Order No.
------------	------	------------------	-------------------------	--------------	--------------

HP-TC BALL / ORAL SURGERY

	HP.TC 2SL	010	1.0	44.5	HP.TC02SL
	HP.TC 4SL	014	1.4	44.5	HP.TC04SL
	HP.TC 6SL	018	1.8	44.5	HP.TC06SL
	HP.TC 8SL	023	2.3	44.5	HP.TC08SL

HP-TC STRAIGHT FISSURE Cross Cut / SURGICAL LENGTH

	HP.TC 557SL	010	3.8	44.5	HP.TC0557SL
	HP.TC 563L	023	5.9	44.5	HP.TC0563L

HP-TC TAPERED FISSURE Cross Cut / SURGICAL LENGTH

	HP.TC 701SL	012	3.8	44.5	HP.TC0701SL
	HP.TC 702SL	016	4.3	44.5	HP.TC0702SL
	HP.TC 703	021	4.8	44.5	HP.TC0703

HP-TC ROUND FISSURE Cross Cut / OPERATIVE

	HP.TC 1557	010	3.8	44.5	HP.TC1557
--	------------	-----	-----	------	-----------

HP-TC DOMED TAPERED FISSURE Cross Cut / ORAL SURGERY

	HP.TC 1702	016	4.4	44.5	HP.TC1702
	HP.TC 1703	021	5.0	44.5	HP.TC1703
	HP.TC 1703L	021	7.0	44.5	HP.TC1703L



D.E.M
Dental Equipment Manufacturer Ltd

Michael House,
Castle Street, Exeter, EX4 3LQ, UK

Phone: +44 (0) 7543 213 432
dem@demdental.com
www.demdental.com



Dental Equipment
Manufacturer Ltd



Carbide Dental Burs

Supporting Excellence in Dentistry



CB-DEM_190213_Lan044-LR