

VirusVisor™

FACE PROTECTION VISOR




COVID-19
PROTECTION



Dental Equipment Manufacturer Ltd



ROBUST PROTECTION WITH THE **VirusVisor™** FACE PROTECTION VISOR



**COVID-19
PROTECTION**



**FULL FACE
PROTECTION**



LIGHTWEIGHT

COMFORTABLE

Soft, protective EPDM lining for sweat reduction

EASY TO ADJUST

Flexible straps secure in any position

WIDE FIELD OF VIEW

With minimal restriction of head movement

NO VAPOUR

Design allows peripheral airflow to remove vapour


EASY CLEANING, REUSABLE

Clean with soap and hot water (50°C), wipe with a soft cloth

LIGHTWEIGHT

Robust yet light PET protective shield, 0.5 mm thick

FITTING



1
Insert the serrated end into the soft band



2
Adjust the head strap to fit comfortably



Dental Equipment Manufacturer Ltd

DEM is a British fully certified medical equipment manufacturer exporting worldwide

DEM Dental Equipment Manufacturer Ltd,
Michael House, Castle Street,
Exeter UK EX4 3LQ
Tel: +44 (0)7746 167589
demdental.com | dem@demdental.com

VirusVisor™

FACE PROTECTION VISOR



COVID-19
PROTECTION



Dental Equipment Manufacturer Ltd

